


Imprint® DNA Modification Kit SIGMA



Company	Input DNA	Time Req'd (hours)
Imprint	>50 pg	<2
Supplier A	>100 pg	>6
Supplier B	>1 µg	>24
Supplier C	>1 ng	>6
Supplier D	>500 ng	>10

Codice Articolo: MOD50-1KT

Breve descrizione del prodotto:

Imprint DNA Modification Kit, For bisulfite DNA conversion & purification

Informazioni aggiuntive:

- Purity Grade: For bisulfite DNA conversion & purification
- URL Link: <https://www.sigmaaldrich.com/US/en/product/SIGMA/MOD50?context=product>
- Image
URL: https://www.sigmaaldrich.com/content/dam/sigma-aldrich/product3/180/mod50_1.eps/_jcr_content/renditions/mod50_1-medium.jpg
- Shipped on ICE/no ICE: STD
- SAP Pack with Pack Qty: 1X1KIT
- Storage Temperature: 20-25C
- Prodotto padre: MOD50-1KT