

REExtract-N-Amp PCR ReadyMix SIGMA



Codice Articolo: R4775-1.2ML

Breve descrizione del prodotto:

REExtract-N-Amp™ PCR ReadyMix™, Ready-to-use 2X PCR Master Mix with Loading Dye

Informazioni aggiuntive:

- Purity Grade: Ready-to-use 2X PCR Master Mix with Loading Dye
- URL Link: <https://www.sigmaaldrich.com/US/en/product/SIGMA/R4775?context=product>
- Image
URL: https://www.sigmaaldrich.com/content/dam/sigma-aldrich/product6/083/r4775_1-2ml_bag.tif/_jcr_content/renditions/r4775_1-2ml_bag-medium.jpg
- Shipped on ICE/no ICE: WICE
- SAP Pack with Pack Qty: 1X1.2ML
- Storage Temperature: -20C
- Prodotto padre: R4775-1.2ML