

## REExtract-N-Amp PCR ReadyMix SIGMA



**Codice Articolo:** R4775-1.2ML

### Breve descrizione del prodotto:

REExtract-N-Amp PCR ReadyMix Ready-to-use 2X PCR Master Mix with Loading Dye

### Informazioni aggiuntive:

- Purity Grade: Ready-to-use 2X PCR Master Mix with Loading Dye
- URL Link: <https://www.sigmaaldrich.com/product/sigma/r4775>
- Image  
URL: <https://api.sigmaaldrich.com/deepweb/assets/sigmaaldrich/product/images/401/063/3f47e17c-9534-4e30-ba6a-b1539b420bc5/640/3f47e17c-9534-4e30-ba6a-b1539b420bc5.jpg>
- Shipped on ICE/no ICE: Wet Ice
- Storage Temperature: -20C
- MSDS Link: <https://www.sigmaaldrich.com/sds/sigma/r4775>
- SAP Pack with Pack Qty: 1x1.2 mL