


REExtract-N-Amp Tissue PCR Kit 10 REACTIONS Sigma-Aldrich



Codice Articolo: XNAT-10RXN

Breve descrizione del prodotto:

REExtract-N-Amp™ Tissue PCR Kit, sufficient for 10 reactions, sufficient for 100 reactions, sufficient for 1000 reactions, hotstart, dNTPs included

Informazioni aggiuntive:

- URL Link: <https://www.sigmaaldrich.com/US/en/product/SIGMA/XNAT?context=product>
- Image
URL: https://www.sigmaaldrich.com/content/dam/sigma-aldrich/product8/165/extract_amp_tissue.ep_s/_jcr_content/renditions/extract_amp_tissue-medium.jpg
- Shipped on ICE/no ICE: WICE
- SAP Pack with Pack Qty: 1X10RXN
- Storage Temperature: -20C
- Prodotto padre: XNAT-10RXN