

Test mixture 2 apolar capillary columns Supelco



Codice Articolo: 86501-10X1ML

Breve descrizione del prodotto:

Test mixture 2 for apolar capillary columns according to Grob, suitable for GC, ready-to-use

Informazioni aggiuntive:

- Purity Grade: suitable for GC, ready-to-use
- URL Link: <https://www.sigmaaldrich.com/US/en/product/SIAL/86501?context=product>
- Image
URL: <https://www.sigmaaldrich.com/deepweb/assets/sigmaaldrich/marketing/global/images/icons/vibrant-m-string/vibrant-m-string.png>
- Shipped on ICE/no ICE: STD
- SAP Pack with Pack Qty: 10X1ML
- Storage Temperature: -2-8°C
- Prodotto padre: 86501-10X1ML